

Product name:Protective respirator

Packing quantity:10pcs/box

Product model:HMS-801

Product style: Non-reusable,folded type

Size:15.5*10.5cm approx.

Material: 5 layers strong protection including:

Skin friendly non-woven fabric, polypropylene meltblown fabric,filtering filling,

Inside non-woven fabric 2 layers,earloop ,nose clip,soft foam,and adjustable hook.

Fitting Instruction

1. Open the protective respirator, face to the inside of the mask, and hold the mask on each hand so that the nose clip is at the top.
2. Hold the protective respirator in position over the nose and mouth, and fasten the other end of the Adjustable hook to the other side of the mask.
3. Adjust to a comfortable position and make the mask fit the face.
4. Bend the nose clip to make a tight seal around the nose.
5. Fit check
 - a) To test the fit of the mask, cup both hands over the mask and inhale sharply. If air flow is felt in the nose area, re-adjust/tighten the noseclip.
 - b) If flows is felt around the edges of the mask, re-position the mask harness to achieve a better fit.
6. Careful observances of these instructions is an important step in safe mask use.

Please read this User Information Sheet carefully before using this product.This product complies with the requirements of EU Regulation (EU)2016/425 for Personal Protective Equipment and meets the requirements of European standard EN 149:2001+A1:2009.FFP2 NR



Check prior to use

The respirator must be selected properly for intended application.An individual risk assessment must be evaluated.Check the respirator that it is undamaged with no visible defects.Check that the expiry date has not been reached(see the packaging).Check the protection class(FFP1 NR/FFP2 NR/FFP3 NR) is appropriate for the product used and its concentration.Do not use the mask if a defect is present or the expiry date has been exceeded.

Intend use of this PPE:

These devices are designed to protect against both solid and liquid aerosols.And be used in industries such as Textile,abattoirs,metallurgy,construction,iron and steel industries,hospitals etc.To Protect against dusts,mists and fumes containing calcium carbonate,clay,kaolin,cellulose,cotton,flour,ferrous,metals,vegetal and mineral oils,metal-working fluids.

This product is designed to protective against the risks:

Risk	Standard Clause	Assessment method
------	-----------------	-------------------

Penetration of particle	EN 149:2001+A1:2009, clause 7.9.1 and 7.9.	Total inward leakage test, Penetration of filter material
Risk to catch fire	EN 149:2001+A1:2009, clause 7.11	Flammability test
Inhalation of CO ₂	EN 149:2001+A1:2009, clause 7.12	Dead space test
Uncomfortable to wear	EN 149:2001+A1:2009, clause 7.13 and 7.16	Head ness check, Breathing resistance

Application/Limitations


This respirator is suitable for use in protection against the non-toxic solid and liquid aerosols. Do not use out of the scope of use defined in the warnings. Failure to properly use this product may result in serious health damage or death.

FFP1 NR: Filter Efficiency 80%; Allocated Protection Factor (FPA) is 4; Examples of applications are Handling of stone/rubble/cellulose.

FFP2 NR: Filter Efficiency 94%; Allocated Protection Factor (FPA) is 10; Examples of applications are Sanding of soft wood, composite materials, rust, putty, plaster, plastics/cutting, deburring, grinding, drilling of metal.

FFP3 NR: Filter Efficiency 99%; Allocated Protection Factor (FPA) is 20; Examples of applications are Sanding of hard wood (beech, oak)/treatment of wood using copper, chrome or arsenic based products/impact stripping of paint/sanding of cement.

Explanation of pictograms

Pictogram	Explanation
	End of shelf life
	See information supplied by the manufacturer
	Temperature range of storage conditions
	Maximum relative humidity of storage conditions

Manufacturing date: 2020-06-16

Warnings

1. Failure to follow all instructions and limitations on the use of this product, or failure to achieve proper fit, may result in damage to your health or death.

2. A properly selected respirator is essential to protect your health. Before using this respirator consult a suitably qualified safety professional to determine the suitability of the product for your intended use.
3. This product does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life. Do not use this respirator when the oxygen concentration is less than 19.5%.
4. Do not use when concentrations of contaminants are immediately dangerous to health or life. Do not use this product in an explosive atmosphere.
5. Leave the work area immediately if: a) breathing becomes difficult or b) dizziness or other distress occurs.
6. Facial hair, beards and certain facial characteristics may reduce the effectiveness of this respirator.
7. Never alter or modify this respirator in any way (except as indicated in the instructions).
8. "NR" means this particle filtering half mask shall not be used for more than one shift. No maintenance is necessary.
9. Respirator should be discarded after use or when it becomes damaged or deformed in any way; no longer forms and effective seal to the face; becomes wet or visibly dirty; breathing becomes difficult.
10. The length of time this respirator can be used depends on contaminants present but should not exceed one shift. The respirator should be replaced sooner if breathing becomes difficult.
11. The head straps are made of synthetic rubber, which may cause skin allergies. Remove the mask immediately if irritation occurs.
12. Keep respirators in the display box away from direct sunlight or contaminants until use. Ambient storage conditions as temperature between -30 °C, to +70°C, and relative humidity <80%.
13. Unless this is fitted according to the "Easy to use" instructions the respirator will not provide the expected level of protection.
14. This respirator is suitable for use in protection against the non-toxic solid and liquid aerosols.
15. Failure to achieve proper fit may result in serious health damage or death.
16. See information supplied by the manufacturer.
17. The respirators must be stored and transported in their original package and protected by the storage temperature and humidity as suggested by the manufacturer.

Manufacturer: Shanghai Homes Mask Co., Ltd

Address: 439# Jinglian Road Shanghai 201108

EU declaration of conformity can be accessed at www.homesbedding.com

Detailed content of the instructions and information notice is the responsibility for manufacturer

Made in China Notified body number: 0598

CE 0598

SGS Fimko Oy, Takomotie 8, FI-00380 Helsinki, Finland.