

Product name: Protective respirator

Packing quantity:10pcs/box

Product model: HMS-902

Product style: Non-reusable,cup type

Material: Strong filtering meltblown fabric,2 layers filtering filling,headband,nose clip,soft foam and valve

Fitting Instruction

1. Cup the respirator in your hand,with the noseclip at your fingertips,allowing the headbands to hang freely below your hand.
2. Position the respirator under your chin with the noseclip up,Pull the top strap over your head at the top back of your head,pull the bottom strap over your head and position it below the ears.
3. Place your fingertips from both hands at the top of the metal noseclip **.Using two hands**,mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of nose piece.
4. Perform a User Seal Check prior to each wearing.
 - a) To test the fit of the respirator, cup both hands over the respirator and inhale sharply. If air flow is felt in the nose area, re-adjust/tighten the noseclip.
 - b) If flows is felt around the edges of the mask, re-position the mask harness to achieve a better fit.
5. Careful observances of these instructions is an important step in safe respirator use.

Please read this User Information Sheet carefully before using this product.This product complies with the requirements of EU Regulation (EU)2016/425 for Personal Protective Equipment and meets the requirements of European standard EN 149:2001+A1:2009.FFP2 NR

Check prior to use



The respirator must be selected properly for intended application.An individual risk assessment must be evaluated.Check the respirator that it is undamaged with no visible defects.Check that the expiry date has not been reached(see the packaging).Check the protection class(FFP1 NR/FFP2 NR/FFP3 NR) is appropriate for the product used and its concentration.Do not use the mask if a defect is present or the expiry date has been exceeded.

Intend use of this PPE:

These devices are designed to protect against both solid and liquid aerosols.And be used in industries such as Textile,abattoirs,metallurgy,construction,iron and steel industries,hospitals etc.To Protect against dusts,mists and fumes containing calcium carbonate,clay,kaolin,cellulose,cotton,flour,ferrous,metals,vegetal and mineral oils,metal-working fluids

This product is designed to protective against the risks:

Risk	Standard Clause	Assessment method
Penetration of particle	EN 149:2001+A1:2009,	Total inward leakage

	clause 7.9.1 and 7.9.	test, Penetration of filter material
--	-----------------------	--------------------------------------

Application/Limitations





This respirator is suitable for use in protection against the non-toxic solid and liquid aerosols. Do not use out of the scope of use defined in the warnings. Failure to properly use this product may result in serious health damage or death.

FFP1 NR: Filter Efficiency 80%; Allocated Protection Factor (FPA) is 4; Examples of applications are Handling of stone/rubble/cellulose.


FFP2 NR: Filter Efficiency 94%; Allocated Protection Factor (FPA) is 10; Examples of applications are Sanding of soft wood, composite materials, rust, putty, plaster, plastics/cutting, deburring, grinding, drilling of metal

FFP3 NR: Filter Efficiency 99%; Allocated Protection Factor (FPA) is 20; Examples of applications are Sanding of hard wood (beech, oak)/treatment of wood using copper, chrome or arsenic based products/impact stripping of paint/sanding of cement.

Explanation of pictograms

Pictogram	Explanation
	End of shelf life
	See information supplied by the manufacturer
	Temperature range of storage conditions
	Maximum relative humidity of storage conditions

Explanation of product marking:

Marking	Description on label	Explanation
	HOMES	Manufacture name, trademark
	HMS-902	Type-identifying marking, style No
	EN 149:2001+A1:2009	Number of European standards
	FFP2 NR	Classification FFP2, non reusable
	CE 0598	CE marking, 7 mm height.

Manufacturing date: 2020-07

Warnings

1. Failure to follow all instructions and limitations on the use of this product, or failure to achieve proper fit, may result in damage to your health or death.
2. A properly selected respirator is essential to protect your health. Before using this respirator consult a suitably qualified safety professional to determine the suitability of the product for your intended use.
3. This product does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life. Do not use this respirator when the oxygen concentration is less than 19.5%.
4. Do not use when concentrations of contaminants are immediately dangerous to health or life. Do not use this product in an explosive atmosphere.
5. Leave the work area immediately if: a) breathing becomes difficult or b) dizziness or other distress occurs.
6. Facial hair, beards and certain facial characteristics may reduce the effectiveness of this respirator.
7. Never alter or modify this respirator in any way (except as indicated in the instructions).
8. "NR" means this particle filtering half mask shall not be used for more than one shift. No maintenance is necessary.
9. Respirator should be discarded after use or when it becomes damaged or deformed in any way; no longer forms and effective seal to the face; becomes wet or visibly dirty; breathing becomes difficult.
10. The length of time this respirator can be used depends on contaminants present but should not exceed one shift. The respirator should be replaced sooner if breathing becomes difficult.
11. The head straps are made of synthetic rubber, which may cause skin allergies. Remove the mask immediately if irritation occurs.
12. Keep respirators in the display box away from direct sunlight or contaminants until use. Ambient storage conditions as temperature between -30 °C, to +70 °C, and relative humidity <80%.
13. Unless this is fitted according to the "Easy to use" instructions the respirator will not provide the expected level of protection.
14. This respirator is suitable for use in protection against the non-toxic solid and liquid aerosols.
15. Failure to achieve proper fit may result in serious health damage or death.
16. See information supplied by the manufacturer.
17. The respirators must be stored and transported in their original package and protected by the storage temperature and humidity as suggested by the manufacturer.

Manufacturer: Shanghai Homes Mask Co., Ltd

Address: 439# Jinglian Road Shanghai 201108

EU declaration of conformity can be accessed at www.homesbedding.com

Detailed content of the instructions and information notice is the responsibility for

manufacturer Made in China Notified body number: 0598  0598

SGS Fimko Oy, Takomotie 8, FI-00380 Helsinki, Finland.